



VET PRODUCTLIST 2024

CRYO IQ



Rev. 2024-01A

The innovative way of Cryosurgery



Phone ▶ +46 31-400 500
Fax ▶ +46 300-56 94 99

E-mail ▶ order@cryoiq.com
Website ▶ www.cryoiq.com



How to order

If possible, please send your order via email.
For customers within the EC community it is important that the VAT number is included in the order.

Order by e-mail: order@cryoiq.com
Order by phone: +46 31-400 500
Order by fax: +46 300-56 94 99

Product information

To find out more about our products, please visit our website www.cryoiq.com.

As a CryoIQ Dealer, we will also provide you with login details to our picture database. Here you will find pictures, artworks, training material, pricelist, etc. all you need for your marketing of CryoIQ products. Please contact our sales department for your login details.

Product number configurations

Learn about product number configuration to make it easier to understand.

CIQ - AA - BB / CC - VET

CIQ CryoIQ Prefix
A Device model
- Divider between device and accessories
B Main Tip / Probe
/ Divider between extra Tip / Probe
C Extra Tip / Probe
VET Veterinary version

Warranty

There is a 1+1 year warranty against technical defects. The warranty is strictly limited to the exchange of defective parts. Dropped devices or broken tips are not under warranty. Other claims for damages such as loss of working hours, incorrect treatment and their consequences, non-executed after treatment and their consequences as well as the non-observance of the safety instructions are excluded from all warranty and liability.

Warranty claims:

All Warranty claims must follow the "Warranty Claim" form that can be downloaded from <http://www.cryoiq.com/warranty>. Please notice that the end-user invoice has to be attached to the Warranty Claim.

Manufacturer



CryoIQ global AB
Apelrödsvägen 1
439 32 Onsala
Sweden

Phone: +46 31-400 500
Fax: +46 300-56 94 99
E-mail: info@cryoiq.com
Website: www.cryoiq.com

Reseller



Devices and options

CryoIQ devices are divided in 2 groups: DERM Plus, PRO depending on functionality and how they could be combined with the various accessories.

DERM plus systems

Order no.	Product description	Pack. size
DERM plus - Devices with fixed tip and gas cartridge with valve		
CIQ-DP-L-VET	CryoIQ DERM Plus Liquid system with 1 device for spray freezing, Ø1-6 mm 1 gas cartridge N ₂ O 25 gram 1 freeze limiter 1 case for storage 1 user guide A6 size	1
BEST SELLER		
CIQ-DP-C1-VET	CryoIQ DERM Plus Contact C1 system with 1 device for contact freezing Ø1 mm 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1
CIQ-DP-C3-VET	CryoIQ DERM Plus Contact C3 system with 1 device for contact freezing Ø3 mm 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1
CIQ-DP-C5 VET	CryoIQ DERM Plus Contact C5 system with 1 device for contact freezing Ø5 mm 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1
CIQ-DP-L/C3/C5-VET	CryoIQ DERM Plus Multi system with 1 device for spray freezing, Ø1-6 mm 1 device for contact freezing Ø3 mm 1 device for contact freezing Ø5 mm 1 freeze limiter 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1

DERM Plus - Spray freezing



DPST

DERM Plus - Contact Freezing

Ø1 mm

Ø3 mm

Ø5 mm



DPC1



DPC3



DPC5



Accessories for DERM Plus devices

Order no.	Product description	Pack. size	Picture
Replacement devices			
CIQ-DAP-DST	Device for spray freezing, Ø1-6 mm	1	DPST
CIQ-DAP-T/C1	Device for contact freezing Ø1 mm	1	DPC1
CIQ-DAP-T/C3	Device for contact freezing Ø3 mm	1	DPC3
CIQ-DAP-T/C5	Device for contact freezing Ø5 mm	1	DPC5
CIQ-DAP-T/C1-5	Contact SET, Ø1, Ø3, Ø5 mm	1	
Accessories			
CIQ-DA-C -VET	Case for storage DERM plus & gas cartridges	1	
CIQ-A-LED	LED light, for 25 gram gas cartridge	1	LED

Gas cartridge 25 gram




25 gram gas cartridge for CryoIQ for devices DERM, DERM Plus and PRO.

Packed in white case



User's Guide A6 size



User's Guide 12 pages, with easy guidelines how to perform treatments and information on clinical studies.



Phone ▶ +46 31-400 500
 Fax ▶ +46 300-56 94 99

E-mail ▶ order@cryoiq.com
 Website ▶ www.cryoiq.com



PRO systems

PRO - Liquid Spray Freezing




PRO - Contact Freezing



User's Guide A6 size

User's Guide 12 pages, with easy guidelines how to perform treatments and information on clinical studies.



Order no.	Product description	Pack. size
PRO - Devices with exchangeable tip and gas cartridge with valve		
CIQ-P-L-VET	CryoIQ PRO liquid system with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 PRO body with lever 1 gas cartridge N ₂ O 25 gram 1 freeze limiter 1 Case for storage 1 user guide A6 size	1
CIQ-P-L/D2-VET	BEST SELLER CryoIQ PRO Liquid D2 system, with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 dermatology tip D2, L=22 mm, spray freezing Ø7-18 mm 1 PRO body with lever 1 gas cartridge N ₂ O 25 gram 1 freeze limiter 1 case for storage 1 user guide A6 size	1
CIQ-P-L/D2/D3-VET	BEST SELLER <i>Used on training video</i> CryoIQ PRO Liquid D2/D3 system with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 dermatology tip D2, L=22 mm, spray freezing Ø7-18 mm 1 dermatology tip D3, L=22 mm, spray freezing Ø9-22 mm 1 PRO body with lever 1 gas cartridge N ₂ O 25 gram 1 freeze limiter 1 case for storage 1 user guide A6 size	1
CIQ-P-C5-VET	CryoIQ PRO Contact C5 system with 1 contact tip for contact freezing Ø5 mm 1 CRYO body with lever 1 gas cartridge N ₂ O 25 gram 1 freeze limiter 1 case for storage 1 user guide A6 size	1
CIQ-P-L/C5-VET	CryoIQ PRO Mix system with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 contact tip for contact freezing Ø5 mm 1 PRO body with lever 1 gas cartridge N ₂ O 25 gram 1 freeze limiter 1 case for storage 1 user guide A6 size	1
CIQ-P-L/D2/C5-VET	CryoIQ PRO DermMix system with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 dermatology tip D2, L=22 mm, spray freezing Ø7-18 mm 1 contact tip for contact freezing Ø5 mm 1 PRO body with lever 1 freeze limiter 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1



CIQ-P-L-VET
 CryoIQ PRO liquid with
 1 cartridge 25 gr




CIQ-P-L/D2/C5-VET
 CryoIQ PRO liquid with Tip
 D2 and C5, 1 cartridge 25 gr




Accessories for PRO devices

Order no.	Product description	Pack. size	Picture
PRO - Accessories			
General			
CIQ-PA-DST60	Standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm	1	PDST
CIQ-PA-DST80	Standard tip 80 µm, L=10 mm, spray freezing Ø7-18 mm	1	PDST
CIQ-PA-BY	PRO body with lever	1	BODY
CIQ-PA-KL	Case for storage PRO & gas cartridge	1	
CIQ-PA-CAP	Protective Cap for PRO glass tip	1	PCAP
CIQ-PA-LMT	Freeze limiter	1	LMT
Dermatology (spray freezing)			
CIQ-PA-D1	Dermatology tip 60 µm, L=22 mm, spray freezing Ø1-6 mm	1	D1
CIQ-PA-D2	Dermatology tip 80 µm, L=22 mm, spray freezing Ø7-18 mm	1	D2
CIQ-PA-D3	Dermatology tip 100 µm, L=22 mm, spray freezing Ø9-22 mm	1	D3
CIQ-PA-D1-3	Dermatology set PRO with 3 tips (D1, D2, D3)	1	
Contact freezing			
CIQ-PA-C1	Contact tip for contact freezing Ø1 mm	1	C1
CIQ-PA-C3	Contact tip for contact freezing Ø3 mm	1	C3
CIQ-PA-C5	Contact tip for contact freezing Ø5 mm	1	C5
CIQ-PA-C7	Contact tip for contact freezing Ø7 mm	1	C7
CIQ-PA-C1-5	Contact set PRO with 3 tips (C1, C3, C5)	1	
Special tip			
CIQ-PA-A45	Angled tip 45°, L=8 cm, spray freezing Ø7-18 mm	1	A45
CIQ-PA-L13	Long tip 100 µm, L=130 mm, spray freezing Ø9-22 mm	1	L13
Other			
CIQ-A-LED	LED light, for 25 gram gas cartridge	1	LED
CIQ-SCALE-A	Portable Scale 500 gram, 0.1 gram	1	SC



LED light with a flexible holder that is EASY to Snap-On our 25 gram gas cartridges.



LMT - Freeze limiter



SC



BODY



Phone ▶ +46 31-400 500
 Fax ▶ +46 300-56 94 99

E-mail ▶ order@cryoiq.com
 Website ▶ www.cryoiq.com



Gas cartridges with OxiBlock® technology

All our gas cartridges are equipped with a valve. The 25 gram gas cartridge fit CryoIQ® models DERM, DERM Plus and PRO, and Cryoalfa® model SUPER, LUX (= CryOmega® Vet). The gas cartridges have a built-in filter 5 µm and are filled with N₂O gas.

Order no.	Product description	Pack. size	Min order qty	Picture
Gas cartridges				
CIQ-G-VN25-VET	Gas cartridge 25 gram N ₂ O with valve for model CryoIQ® DERM, DERM Plus, PRO and Cryoalfa® SUPER, LUX (version 2).	1	10 single package	G4
CIQ-G-VN25x4-VET	New! Gas cartridge 25 gram N ₂ O with valve for model CryoIQ® DERM, DERM Plus, PRO and Cryoalfa® SUPER, LUX (version 2)	4	5 (4-pack)	
CIQ-G-VN25x10-VET	New! Gas cartridge 25 gram N ₂ O with valve for model CryoIQ® DERM, DERM Plus, PRO and Cryoalfa® SUPER, LUX (version 2)	10	1 (10-pack)	
CIQ-G-VN25-LUX-VET	Gas cartridge 25 gram N ₂ O with valve for model Cryoalfa® LUX (version 1 and 2) and CryOmega® Vet	1	10 single package	

*Cryoalfa LUX (Version 2) in general produced 2012 or later, most of LUX version 2 has silver lever.



G1 / With valve
25 gram version



G4 / With valve
25 gram version for
Cryoalfa LUX



1-pack (G1)



4-pack (G2)



10-pack (G3)

OxiBlock®

The new OxiBlock® coating on the inside of the cartridge gives you a contamination-free environment for the gas. What we have seen is that without our OxiBlock® coating you are not able to obtain a contamination inside of the gas cartridge. This affects the durability and purity of the gas.

Warranty PLUS

With the meticulous control of all production steps of the gas cartridges, we are able to offer you, that if any faults in the devices from CryoIQ® or Cryoalfa® appears caused by use of CryoIQ® gas cartridges we will replace unit with a new equal CryoIQ® device. Just remember that it is important that you follow the cleaning instructions in our User Guide.

Filling

All of our cartridges are filled in-house. We have put a lot of efforts in developing a unique gas filling station that has the latest technique to optimize filling of cartridges. To establish a more precise and clean filling process, the filling station itself is equipped with nine separate 5 micron gas filters to ensure 100% purity. With our in-house filling process, we have 100% traceability of the gas used in our cartridges.

It was important for us to find an independent solution where we could fill and refill our cartridges. We have a close cooperation with Nippon Gas who is our consultants and manufacturer of gas, but more importantly, we have the in-house knowledge. Today we have in our production two filling stations and are now looking for a fully automated filling line, with a capacity of more than 500 gas cartridges per day to meet future demand.

N₂O gas

With our filling process, we can ensure that the gas we use is kept 100% clean during the entire filling process.

Cap

The fine and precise thread in the valve requires a precise cap to seal completely. Our new cap construction is in fine metal and it has its own sealing to ensure an even better shelf life for the cartridges.

Shelf Life

Together with our gas manufacturer, Praxair, we have implemented a Just-In-Time system so that we can always offer you as our dealer the best shelf life for your gas cartridges. From the production day of the gas the durability is 3 years and with our system we never have gas older than 1 month in our production.

What we have seen before is that larger batches of production results in shorter shelf time for the gas cartridges. With our Just-In-Time system we have adjusted our production so that we never need to hold unnecessary stock. This system is also very cost effective which enables us to keep our prices low.




Product pictures

To make life easy for you, all our pictures and leaflets are available for download via our picture database. <http://www.picture.cryoiq.com>

Send us an e-mail to info@cryoiq.com and we will create a personal login for you.

If you need any specific picture or would like to have designed picture. Ask us and we will see how we can support you.



My notes




Phone ▶ +46 31-400 500
 Fax ▶ +46 300-56 94 99

E-mail ▶ order@cryoiq.com
 Website ▶ www.cryoiq.com




Marketing material

We offer following marketing material in smaller volume, for large quantities please contact us.
 In our graphic studio we can produce custom designed material.



Veterinary

Description: Leaflets for CryoIQ products
 Format: Printed
 Languages: English, Russian
 How to order: CIQ-F-MED-VET-EN (RU)
 Option: Can be printed with Partner logo
 Cost: Depending on volume



Userguide
 veterinary

Description: Userguides for CryoIQ products
 Format: Printed
 How to order: B-UG-MED-EN
 Languages: Bulgarian, English, Italian, Russian, Swedish,
 more languages on the way.



Description: Rollup showing product portfolio
 Format: RollUp 80 x 200 cm
 How to order: CIQ-PRINT-ROLLPORT2016
 Languages: Universal
 Cost: €80 (approximately)



Description: DVD with pictures and other marketing material
 Format: DVD
 How to order: CIQ-DVD-2024
 Languages: -
 Cost: -